

40.401

HYGIENE PP LAB COAT 40 g/m²

White with two pockets and snap fastener

FEATURES & BENEFITS

- Single use hygiene lab coats protect the working environment from particles coming from the body of the operator (dander, hair, sweating ...)
- They allow the operator to protect himself from the risks of soiling related to the activity he carries out.

PRODUCT DETAILS

Product name:	Hygiene PP lab coat SafeWear
Product type :	Single use, non sterile
Material :	Polypropylene
Confection :	Sewed & elasticated (collar and wrists) Snap fastener
Color :	White
Grammage :	40 g/m ² (± 10%)
Packaging :	50 individual polybags
Origin :	China



Item	Size	Length (± 2 cm)	Width (± 2 cm)
40.401/02	M	104 cm	130 cm
40.401/03	L	104 cm	135 cm
40.401/04	XL	104 cm	140 cm
40.401/05	XXL	104 cm	145 cm
40.401/06	XXXL	104 cm	150 cm

CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001.

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICS

OUTERCASE SPECIFICATIONS					PACKAGING SPECIFICATIONS	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	PACKAGING	EAN CODE
40.401/02		4,1		3 700 381 505 596		3 700 381 505 602
40.401/03		4,6		3 700 381 500 829		3 700 381 502 274
40.401/04	470 x 280 x 280	4,9	36 cases	3 700 381 500 836	50 Individual polybags	3 700 381 502 281
40.401/05		5,3		3 700 381 500 843		3 700 381 502 298
40.401/06		5,8		3 700 381 506 722		3 700 381 506 739