

J50656

CPE APRON GOWN SPLASH RESISTANT

FEATURES & BENEFITS

- Provides superior protection to a traditional apron
- Covers almost the whole body but also the arms
- Resistant PE tie bonds.
- « Thumb loops » system to avoid the sleeves go back
- It adapts to different morphologies

TECHNICAL DESCRIPTION

Product name : CPE apron gown splash resistant
 Product type : Single use
 Material : CPE
 Confection : Cutting, folding, sewing
 Color : Light blue
 Thickness : 35 µm
 Dimensions : 119 x 105 cm
 Unit weight : 54 g
 Manufacturing country : China



CERTIFICATIONS & STANDARDS

Manufacturing site : ISO 9001.

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICS

OUTER CASE SPECIFICATIONS					PACKAGING SPECIFICATIONS	
Article	SIZE (mm)	WEIGHT (kg)	Qty / pallet	EAN CODE	Packaging	EAN CODE
J50656	380 x 200 x 310	3,9	48 cases	3 662 036 019 961	3 bags of 25 pcs	3 662 036 019 707