

40.298 SAFEWEAR SMS ROBE

40.299 Long sleeves

FEATURES & BENEFITS

- **Opaque SMS fabric** : Lightweight and durable material providing comfort while preserving patient privacy.
- **Loose fit with waist tie** : Adjustable fit for secure coverage and easy dressing before or after a procedure.
- **Wearing comfort**: Soft and breathable fabric suitable for hospital environments.
- **Practical design** : Simple opening for easy dressing and quick access for care.



PRODUCT DETAILS

Product name :	SafeWear SMS Robe
Material :	SMS
Grammage :	35 g/m ²
Confection :	Sewn Long sleeves
Fastener :	Waist tie
Seams :	Overlocked
Qty / case :	100 units
Origin :	China

Item	Color	Size (cm)
40.298	Blue 	110 x 145
40.299	Blue 	125 x 150

CERTIFICATIONS & STANDARDS

Manufactured under an ISO 9001 certified quality management system.

STORAGE CONDITIONS

Normal storage and conservation conditions: Keep away from moisture and direct sunlight. Store at room temperature.

40.298 SAFEWEAR SMS ROBE

40.299 Long sleeves

LOGISTICS



Outer case
spec.

Item	Size Mm	QTY /pallet	Packaging
40.298	550 x 270 x 290	20 cases (5 layers of 4 cases)	10 polybags of 10 units
40.299	650 x 280 x 330		



EAN Codes

Item	Outer case	Polybag
40.298	6973122514288	6973122512987
40.299	6973122514295	6973122512994