

40.461

SAFEWEAR ADULT PATIENT GOWN

Without pockets – Short sleeves
Medical Device Class I

FEATURES & BENEFITS

- **Soft and opaque polypropylene fabric:** Antistatic and skin-friendly material that preserves privacy while helping to reduce the risk of contamination.
- **Loose fit designed for adults :** Makes dressing easier and ensures greater comfort for patients before surgery.
- **Secure fastening:** Neck and waist ties for easy and reassuring adjustment without restricting movement.
- **Hygienic use :** Single-use design helping to reduce the risk of healthcare-associated infections in hospital environments.

PRODUCT DETAILS

Product name :	SafeWear adult patient gown
Material :	Spunbond polypropylene (PP)
Grammage :	42 g/m ²
Confection :	Sewn
	Short sleeves
	Round neck
Fastener :	Back opening with neck and waist ties
Seams :	Overlocked
Qty / case :	100 units
Origin :	China



Item	Color	Size (cm)
40.461/02	Bleu 	110 x 135
40.461/03	Bleu 	125 x 140

Données de résistance mécanique

Test	Standard	Laboratory	Report N°	Result
Grammage	NF EN ISO 9073-1	Cetelor	25-0321-1	42 g/m ² (+/- 5%)
Tensile strength	NF EN ISO 9073-3	Cetelor	25-0321-1	Machine direction: 73.7 N Cross direction: 55.3 N
Elongation	NF EN ISO 9073-3	Cetelor	25-0321-1	Machine direction: 68.06% Cross direction: 74.98%

CERTIFICATIONS & STANDARDS

- Compliant with the requirements of Regulation (EU) 2017/745 on medical devices.
- Manufactured under an ISO 9001 & ISO 13485 certified quality management system

STORAGE AND CONSERVATION CONDITIONS

Normal storage and conservation conditions: Keep away from moisture and direct sunlight. Store at room temperature.

LOGISTICS

Outer case
spec

Item	Size mm	QTY / pallet	Packaging
40.461/02	460 x 290 x 440	18 cases (3 layers of 6 cases)	10 polybags of 10 units
40.461/03	480 x 300 x 460		



EAN Codes

Item	Outer case	Polybag
40.461/02	6973122517036	6973122517029
40.461/03	6973122514035	6973122516138