

40.803 SCRUB SUIT : TUNIC & PANTS

40.804

With 3 pockets

40.805

Medical Device Class I

FEATURES & BENEFITS

- **SMS material:** Three-layer construction providing an excellent balance of mechanical strength, lightweight comfort and breathability.
- **Optimal protection:** Effective barrier against particles, fluids and contaminants, helping to reduce the risk of cross-contamination in the operating room.
- **Long-lasting comfort:** Soft, non-irritating, antistatic and breathable fabric, allowing healthcare professionals to work comfortably for extended periods.
- **Freedom of movement:** Ergonomic design providing excellent freedom of movement during technical procedures and everyday activities.
- **Practical storage:** Convenient storage space for essential accessories used in the care environment.



PRODUCT DETAILS

Product name :	SafeWear scrub suit
Material :	SMS
Grammage :	35 g/m ²
Confection :	Sewn
	Shorts sleeves
	V-neck with binding
	With 3 pockets
	Waist tie
	Elasticated ankles
Seams :	Overlocked
Qty / case :	50 units – individual packaging
Origin :	Myanmar

Item	Color
40.803	 Blue
40.804	 Green
40.805	 Mauve

CERTIFICATIONS & STANDARDS

- Compliant with the requirements of Regulation (EU) 2017/745 on medical devices.
- Manufactured under an ISO 9001 & ISO 13485 certified quality management system.

STORAGE CONDITIONS

Normal storage and conservation conditions: Keep away from moisture and direct sunlight. Store at room temperature.

40.803
40.804
40.805

SCRUB SUIT : TUNIC & PANTS

With 3 pockets
Medical Device Class I

Mechanical resistance data				
Test	Standard	Laboratory	Report n°	Result
Grammage	NF EN ISO 9073-1:2023	Cetelor	25-0359-4	35 g/m ² (+/- 6%)
Tensile strength	NF EN ISO 9073-3:2023	Cetelor	25-0359-4	Machine direction: 79.1 N Cross direction: 40.5 N
Elongation	NF EN ISO 9073-3	Cetelor	25-0359-4	Machine direction: 74.2% Cross direction: 68.1%
Thickness	NF EN ISO 9073-2:1997	Cetelor	25-0359-4	0,21 mm
Thermal resistance	NF EN ISO 11092:2014	Cetelor	25-0359-4	0,019 m ² .K/W
Water vapour resistance	NF EN ISO 11092:2014	Cetelor	25-0359-4	1,15 m ² .Pa/W

NF EN 13795-2: Surgical clothing and drapes – Requirements and test methods – Part 2: Surgical clothing

Tests	Laboratory	Report n°	Required level	Result
Bursting strength (dry condition)	SGS	SL52405251351301TX	≥ 40 kPa	145 kPa
Tensile strength (dry condition)			≥ 20 N	Machine direction: 48.6 N Cross direction: 89.6 N
Particle release			≤ 4,0	Face: 2.19 Back: 2.27
Microbial penetration (dry condition)			≤ 100 CFU	Compliant

40.803
40.804
40.805

SCRUB SUIT : TUNIC & PANTS

With 3 pockets
Medical Device Class I

LOGISTICS



Outer case spec.

Item	Product size	Case size	QTY/ Pallet	Packaging
40.803/XS	XS	460 x 300 x 280	30 cases (5 layers of 6 cases)	50 units / case
40.803/01	S	460 x 330 x 300	25 cases (5 layers of 5 cases)	
40.803/02	M	460 x 330 x 300		
40.803/03	L	460 x 330 x 300		
40.803/04	XL	480 x 330 x 370	16 cases (4 layers of 4 cases)	
40.803/05	XXL	480 x 330 x 370		
40.803/06	XXXL	480 x 330 x 370		
40.804/XS	XS	460 x 300 x 280	30 cases (5 layers of 6 cases)	50 units / case
40.804/01	S	460 x 330 x 300	25 cases (5 layers of 5 cases)	
40.804/02	M	460 x 330 x 300		
40.804/03	L	460 x 330 x 300		
40.804/04	XL	480 x 330 x 370	16 cases (4 layers of 4 cases)	
40.804/05	XXL	480 x 330 x 370		
40.804/06	XXXL	480 x 330 x 370		
40.805/XS	XS	460 x 300 x 280	30 cases (5 layers of 6 cases)	50 units / case
40.805/01	S	460 x 330 x 300	25 cases (5 layers of 5 cases)	
40.805/02	M	460 x 330 x 300		
40.805/03	L	460 x 330 x 300		
40.805/04	XL	480 x 330 x 370	16 cases (4 layers of 4 cases)	
40.805/05	XXL	480 x 330 x 370		
40.805/06	XXXL	480 x 330 x 370		



EAN Codes

Item	Outer case	Item	Outer case	Item	Outer case
40.803/XS	8834000180378	40.804/XS	8834000180514	40.805/XS	8834000180651
40.803/01	8834000180392	40.804/01	8834000180538	40.805/01	8834000180675
40.803/02	8834000180415	40.804/02	8834000180552	40.805/02	8834000180699
40.803/03	8834000180439	40.804/03	8834000180576	40.805/03	8834000180712
40.803/04	8834000180453	40.804/04	8834000180590	40.805/04	8834000180736
40.803/05	8834000180477	40.804/05	8834000180613	40.805/05	8834000180750
40.803/06	8834000180491	40.804/06	8834000180637	40.805/06	8834000180774