

50.805
50.806

SAFEWEAR WATERPROOF & ABSORBENT APRON

Pulp + PE

Dispenser bag

FEATURES & BENEFITS

- **Dual-layer construction:** absorbent pulp layer to retain liquids and waterproof PE film preventing liquid penetration.
- **Effective protection:** Ideal for procedures involving splashes or fluids, helping to keep clothing dry.
- **Easy to put on:** head-opening design for quick and hassle-free donning.
- **Bag packaging:** Ensures hygienic presentation and helps minimize the risk of contamination before use.



PRODUCT DETAILS

Product name :	SafeWear waterproof & absorbent apron
Material :	Pulp + polyethylene (PE)
Confection :	Cutting / Folding
Grammage :	Pulp 23 g/m ² Polyethylene 17 g/m ²
Packaging :	Dispenser bag
Origin :	China

Item	Color	Size
50.805	Blanc / Bleu	120 x 70 cm
50.806	Blanc / Bleu	135 x 70 cm

CERTIFICATIONS & STANDARDS

Manufactured under an ISO 9001 certified quality management system.

Mechanical Resistance Data

Test	Standard	Laboratory	Report N°	Result
Grammage (Pulp+PE)	NF EN ISO 12625-6:2017	Cetelor	25-0321-3	40 g/m ² (+/- 10%)
Tensile strength	NF EN ISO 12625-4:2017	Cetelor	25-0321-3	Machine direction: 32.8 N Cross direction: 25.9 N
Elongation at break	NF EN ISO 12625-4:2017	Cetelor	25-0321-3	Machine direction: 62.8% Cross direction: 41.4%
Tensile energy absorption at break	NF EN ISO 12625-4:2017	Cetelor	25-0321-3	Machine direction: 122 J/m ² Cross direction: 59 J/m ²

50.805
50.806

SAFEWEAR WATERPROOF & ABSORBENT APRON
Pulp + PE
Dispenser bag

STORAGE CONDITIONS

Normal storage and conservation conditions: Keep away from moisture and direct sunlight. Store in a cool or room-temperature environment.



Outer case spec

Item	Size Mm	QTY/ Pallet	Packaging
50.805	450 x 370 x 180	36 cases (9 layers of 4 cases)	4 polybags of 50 pcs
50.806	500 x 370 x 180		

EAN Codes

Item	Outer case	Polybag
50.805	6973122517050	6973122517043
50.806	6973122517074	6973122517067