

70.722
70.732

SAFEWEAR AERATED TIE-BACK SURGICAL HOOD

Medical Device Class I

FEATURES & BENEFITS

- Polyester and viscose non-woven fabric
- Air-permeable material
- Highly covering design
- LATEX-FREE
- Highly breathable for enhanced comfort
- High stretchability – suitable for all hair types and volumes



RECOMMENDED USE

- Suitable for use by healthcare professionals in operating rooms and care units.

PRODUCT DETAILS

Product name :	Safewear aerated tie-back surgical hood
Material :	Polyester & viscose
Grammage :	40g/m ²
Dimensions :	Height 25cm (±2 cm) Length 25cm (±1 cm)
Packaging :	Dispenser Box of 100 pcs
Manufacturing country :	China

Item	Color
70.722	 Blue
70.732	 Green

MECHANICAL INFORMATION

Test	Standard	Laboratory	Report n°	Results
Grammage	NF EN ISO 9073-1: 2023	Internal laboratory	B27250916B07	40 g/m ² ± 5%
Tensile strength	NF EN ISO 9073-3: 2023			MD: 58 N/5cm ±10% CD: 35 N/5cm ±10%
Elongation	NF EN ISO 9073-3: 2023			MD: 36 % ±10% CD: 78 % ±10%

70.722

70.732

SAFEWEAR AERATED TIE-BACK SURGICAL HOOD

Medical Device Class I

CERTIFICATIONS & STANDARDS

Meets the requirements of Regulation (EU) 2017/745 on medical devices.

Manufacturing site : ISO 9001 and ISO 13485 certified.

PRECAUTIONS

The product must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation on the subject of disposal of this product once it has been used.

STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

LOGISTICAL INFORMATION

Outer carton				Inner Box
Item	Size mm	Qty/ pallet	Packaging	Size mm
70.722	470 x 350 x	12 cartons	6 boxes of 100 pcs	340 x 150 x 230
70.732	475	(3 layers of 4 cartons)		

EAN Codes

Item	Carton	Inner box
70.722	6973122517111	6973122517104
70.732	6973122517722	6973122517326