

70.942  
70.943

## SAFEWEAR TIE-BACK SURGICAL CAP

Medical Device Class I

### FEATURES & BENEFITS

- 100% polypropylene non-woven fabric
- Ties to be fastened at the nape of the neck
- LATEX-FREE
- Soft and lightweight material
- Highly breathable for enhanced comfort
- High stretchability – suitable for all hair types and volumes



### RECOMMENDED USE

- Suitable for use by healthcare professionals in operating rooms and care units.

### PRODUCT DETAILS

|                         |                                              |
|-------------------------|----------------------------------------------|
| Product name :          | Safewear PP tie-back surgical cap            |
| Material :              | Polypropylene                                |
| Grammage :              | 25g/m <sup>2</sup>                           |
| Dimensions :            | Height 49 cm (±2 cm)<br>Length 22 cm (±2 cm) |
| Packaging :             | Dispenser Box of 100 pcs                     |
| Manufacturing country : | China                                        |

| Item   | Color                                                                                       |
|--------|---------------------------------------------------------------------------------------------|
| 70.942 |  Blue  |
| 70.943 |  Green |

### MECHANICAL INFORMATION

| Test             | Standard               | Laboratory          | Report n°    | Results                                |
|------------------|------------------------|---------------------|--------------|----------------------------------------|
| Grammage         | NF EN ISO 9073-1: 2023 | Internal laboratory | B19250314B06 | 25 g/m <sup>2</sup> ± 5%               |
| Tensile strength | NF EN ISO 9073-3: 2023 |                     |              | MD: 47 N/5cm ±10%<br>CD: 34 N/5cm ±10% |
| Elongation       | NF EN ISO 9073-3: 2023 |                     |              | MD: 103 % ±10%<br>CD: 107 % ±10%       |

70.942

# SAFEWEAR TIE-BACK SURGICAL CAP

70.943

Medical Device Class I

## CERTIFICATIONS & STANDARDS

Meets the requirements of Regulation (EU) 2017/745 on medical devices.

Manufacturing site : ISO 9001 and ISO 13485 certified.

## PRECAUTIONS

The product must be used on intact skin only.

The reuse or extended use may lead to infection or cross-contamination.

Please refer to your national legislation on the subject of disposal of this product once it has been used.

## STORAGE CONDITIONS

Normal conditions of conservation & storage : not to be exposed to moisture and sunlight, must be stored in a cool environment or at room temperature.

## LOGISTICAL INFORMATION

| Outer carton |                 |                         |                     | Inner Box       |
|--------------|-----------------|-------------------------|---------------------|-----------------|
| Item         | Size mm         | Qty/<br>pallet          | Packaging           | Size mm         |
| 70.942       | 670 x 280 x 415 | 12 cartons              | 10 boxes of 100 pcs | 270 x 200 x 130 |
| 70.943       |                 | (3 layers of 4 cartons) |                     |                 |

## EAN Codes

| Item   | Carton        | Inner box     |
|--------|---------------|---------------|
| 70.942 | 6973122517098 | 6973122517081 |
| 70.943 | 6973122517937 | 6973122519436 |