Medicom

Technical Data Sheet

40.420
40.421

SAFEWEAR WATERPROOF ISOLATION / PATIENT CARE GOWN

Jersey cuffs
Medical Device Class I

FEATURES & BENEFITS

- **PE-coated PP fabric** : Waterproof material providing excellent resistance to liquids and reliable protection for healthcare professionals.
- **Enhanced protection** : Long sleeves and a full-coverage design providing an effective barrier during procedures involving exposure to fluids.
- **Jersey cuffs** : Soft, fitted cuffs allowing gloves to be securely positioned over the gown.
- **Back closure**: Reduces frontal exposure and ensures a secure fit throughout procedures.
- **Optimal hygiene** : Single-use design helping to reduce the risk of cross-contamination.



Item	Color
40.420	White
40.421	Blue

PRODUCT DETAILS

Product name :	**SafeWear waterproof isolation gown**
Material :	Polypropylene (PP) + Polyethylene (PE)
Grammage :	30 g/m²
Confection :	Sewn Round neck Long sleeves Jersey cuffs
Fastener :	Waist and neck ties
Dimensions :	110 x 135 cm
Seams :	Overlocked
Qty/case :	100 units
Origin:	China

Données de résistance mécanique				
Test	**Standard**	**Laboratory**	**Report N°**	**Result**
Grammage	NF EN ISO 9073-1	Cetelor	25-0321-2	30 g/m² (+/- 7%)
Tensile strength	NF EN ISO 9073-3		25-0321-2	**Machine direction:** 41.8 N **Cross direction:** 20.8 N
Elongation	NF EN ISO 9073-3		25-0321-2	**Machine direction:** 62.00% **Cross direction:** 55.86%
Resistance to water penetration	NF EN ISO 811		2602_0013	50 485 Pa

SAFEWEAR WATERPROOF ISOLATION / PATIENT CARE GOWN

Jersey cuffs
Medical Device Class I

CERTIFICATIONS & STANDARDS

Compliant with the requirements of Regulation (EU) 2017/745 on medical devices.
Manufactured under an ISO 9001 & ISO 13485 certified quality management system.

STORAGE AND CONSERVATION CONDITIONS

Normal storage and conservation conditions: Keep away from moisture and direct sunlight. Store at room temperature.

LOGISTICS

Outer case spec.

Item	Size Mm	QTY/ pallet	Packaging
40.420 **40.421**	500 x 330 x 400	12 cases (3 layers of 4 cases)	10 polybags of 10 units

EAN codes

Item	Outer case	Polybag
40.420	6973122514608	6973122514202
40.421	6973122516992	6973122516985

Medicom SAS
Bld de la chanterie | 49124 Saint Barthélemy d'Anjou – France
T : 33 (0) 241 963 434 | F : + 33 (0) 241 963 453 | sales.support.pro.eu@medicom.com | www.medicom-eu.com